

**5X TAPS-PEG 8000 Buffer for Tn5 Tagmentation:**

50 mM TAPS-NaOH at pH 8.5

25 mM MgCl<sub>2</sub>

40% PEG 8000

**Notes:**

Dissolving the PEG 8000 in may take several hours. It is recommended to allow this to happen overnight at room temperature on a stir plate. Shaking is likely insufficient to fully dissolve the PEG 8000.

The TAPS concentration has been reduced from 1.04 M to 50 mM. If planning to use the original concentration, add an additional 24.064 g TAPS to the TAPS-MgCl<sub>2</sub> buffer (25.28 g Total).

**Reagents:**

| Chemical          | Concentration | CAS                  | URL   |
|-------------------|---------------|----------------------|---|
| TAPS              | 243.28 g/mol  | 29915-38-6           | <a href="https://us.vwr.com/store/catalog/static_catalog.jsp?catalog_number=97064-208">https://us.vwr.com/store/catalog/static_catalog.jsp?catalog_number=97064-208</a>   |
| MgCl <sub>2</sub> | 1M            | 7791-18-6            | <a href="https://www.gbiosciences.com/Buffers-Reagents-Chemicals/Molecular-Biology-Related-Buffers-Chemicals/Magnesium-Chloride-1M">https://www.gbiosciences.com/Buffers-Reagents-Chemicals/Molecular-Biology-Related-Buffers-Chemicals/Magnesium-Chloride-1M</a> |
| NaOH              | 10N           | 1310-73-2, 7732-18-5 | <a href="https://www.fishersci.com/shop/products/sodium-hydroxide-solution-10n-certified-fisher-chemical-3/SS2551">https://www.fishersci.com/shop/products/sodium-hydroxide-solution-10n-certified-fisher-chemical-3/SS2551</a>                                   |
| PEG 8000          | 7000 to 9000  | 25322-68-3           | <a href="https://www.fishersci.com/shop/products/polyethylene-glycol-8000-peg-fisher-bioreagents-2/BP233100">https://www.fishersci.com/shop/products/polyethylene-glycol-8000-peg-fisher-bioreagents-2/BP233100</a>   |

**Protocol:****To a 200 ml autoclaved beaker:**

- Makes 50 mls TAPS-MgCl<sub>2</sub> buffer.
1. Add an autoclaved stir-bar.
  2. Add 1.2164 g TAPS.
  3. Add 30 mls Ultra-Pure H<sub>2</sub>O.
  4. Add 2.5 mls 1M MgCl<sub>2</sub>.
  5. Stir until TAPS is fully dissolved.
  6. Adjust to pH 8.5 with 10N NaOH.
  7. Transfer to 50 ml Falcon and adjust final volume of TAPS-MgCl<sub>2</sub> buffer to 50 mls with Ultra-Pure H<sub>2</sub>O.

Filter sterilize using .22 um filter (Optional, recommended if pH probe contamination is a concern).

**To a 200 ml autoclaved beaker:**

- Makes 100 mls 5X TAPS-PEG 8000 using TAPS-MgCl<sub>2</sub> buffer made above.
1. Add an autoclaved stir-bar.
  2. Add 40 grams PEG 8000.
  3. Add 50 mls of TAPS- MgCl<sub>2</sub> solution to 50ml falcon.
  4. Add 10 mls Ultra-Pure H<sub>2</sub>O
  5. Cover and stir until solution is clear and homogeneous. (May take Overnight)
  6. Transfer to falcon tubes.